



Changing the face of sleep

Wisp



PHILIPS

RESPIRONICS

Congratulations

You are about to experience our latest innovation in minimal contact mask technology, Wisp.

Wisp, a minimal contact nasal mask designed around you.

- Wisp provides an open field of vision for feeling of greater independence
- The cushion sits below your nose bridge, making minimal contact with your face while creating an excellent seal
- Multiple cushion sizes included with the mask to achieve your best fit



Multiple cushion sizes included with the mask to achieve your best fit.

Wisp is designed to comfortably fit more than 98% of users.*

*2012 Wisp customer preference trial, data on file



Wisp features just for you:

- Contemporary frame options to complement your style
- Soft and comfortable for an exceptional quality of sleep
- Easy to adjust for an excellent seal and personalized custom fit


Self-adjusting strap with a loop that can hold mask tubing when the tube is placed over the top of your head

Tip-of-the-nose cushion design

Clips make the headgear easy to adjust

Reversible fabric frame; also available in silicone

Elbow quick release tabs



Helpful tips

Here are some helpful tips to keep in mind when using your Wisp minimal contact nasal sleep apnea mask:



- 1 Your Wisp cushion may have a different airflow sensation because it's designed to fit closer to your nose. If you don't feel like you're getting enough air, try one of the larger size cushions included with your mask.
- 2 Do not over-tighten the mask to create a seal. Loosen the straps to allow the cushion to inflate with air and create a comfortable seal on your face.
- 3 Wisp was designed to cover less of your face while not having cushions inserted into your nose. The Wisp cushion fits around the tip of your nose to deliver therapy.

Additional resources

We're dedicated to your satisfaction and aim to answer all of your questions about the Wisp sleep apnea mask. For more information, including helpful videos, visit www.philips.us/yourwisp



Call us

Patient Management Center
1-866-820-3548

Philips Respironics Wisp is available by prescription only. As with any mask or interface, users may experience certain types of discomfort or irritation. Talk to your doctor to determine if it is right for you. Be sure to ask your homecare provider for more information about Wisp.

PHILIPS


RESPIRONICS

**Philips Healthcare is part of
Royal Philips**

How to reach us

www.philips.com/healthcare

Philips Respirationics
1010 Murry Ridge Lane
Murrysville, PA 15668

Respirationics and Wisp are trademarks of Respirationics, Inc. and its affiliates.

Please visit www.philips.us/yourwisp



© 2015 Koninklijke Philips N.V. All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from use of this publication. CAUTION: U.S. federal law restricts these devices to sale by or on the order of a physician.

RRDPGH MG 06/01/15 MCI 4106728 PN 1109130